Tensoval® duo control

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■ English

- 1 Cuff socket
- 2 Connection to computer (USB)
- 3 Mains adapter socket
- 4 Extra large LCD display
- 5 Identification fields for user names
- 6 Memory button for user 2
- 7 Memory button for user 1
- 8 START/STOP button
- 9 Battery compartment with mechanism for easy battery removal
- 10 Grip tab for pulling the cuff on
- 11 Sizing scale for correct adjustment of the cuff
- 12 Cuff connector with extra wide grip
- 13 Extra long tube
- 14 Drawing to help with the correct applying of the cuff
- 15 Recess for bend of the elbow: ergonomic shape of cuff prevents incorrect applying and slipping
- 16 Easy to open Velcro fastener

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Preliminary remarks

Please read these instructions for use carefully before first use, as correct blood pressure measurement depends on the appropriate use of the device.

These instructions for use are designed to instruct you, from the very start, in the individual steps of self-measurement of blood pressure using Tensoval duo control. You will thus receive important and helpful hints for producing reliable results for your personal blood pressure profile. Use this device according to the information described in the instructions for use and keep these in a safe place.

Duo Sensor Technology

The innovative Duo Sensor Technology combines two professional measuring technologies: the oscillometric technology and the Korotkoff technology.

Whilst most automatic blood pressure measurement devices only use oscillometric technology, Duo Sensor Technology also utilises Korotkoff's highly accurate method of blood pressure measurement, which is also used by physicians to measure blood pressure. It is characterised by its reduced susceptibility to malfunction, producing correct measured values even in patients with heart rhythm disorders.

Physicians use a stethoscope to listen to the so-called Korotkoff sounds and thereby to determine blood pressure; the Tensoval duo control does exactly the same via an integrated microphone.

Duo Sensor Technology thus produces accurate results and is still easy to use.

Comfort Air Technology

The Comfort Air Technology enables approximate determination of the systolic blood pressure value during inflation and on that basis determines the individually required inflation pressure for blood pressure measurement. This allows measurements on the upper arm to be taken with less discomfort.

Table of contents

1. Introduction	92
2. General information on blood pressure	92
2.1 Significance of blood pressure values	92
2.2 Importance of self-measurement of blood pressure	92
3. Getting ready for self-measurement	93
3.1 Inserting / changing the batteries	93
3.2 Setting date and time	93
3.3 10 golden rules for blood pressure measurement	94
3.4 Application of the cuff	94
4. Measuring the blood pressure	95
5. Memory function	96
5.1 Storing measured values	96
5.2 Recalling the measured values	96
5.3 Deleting the measured values	97
5.4 Operation of guest mode	98
6. Explanation of error displays	98
7. Explanation of symbols	100
8. Safety information concerning your health	100
8.1 Instructions for self-measurement	100
8.2 Arrhythmias, heart rhythm disorders	101
9. Maintenance of the device	101
10. Accessories and spare parts	102
11. Transmission of blood pressure values to a computer	102
12. Warranty conditions	102
13. Technical data	103
14. Power supply, disposal instructions	104
14.1 Batteries, mains adapter	104
14.2 Disposal	104
15. Safety information concerning the device	104
16. Legal requirements and guidelines	104
17. Instructions for measurement function check	105
18. Instructions for the calibration mode	105
19. Contact information for customer queries	105

Category	Systolic pressure	Diastolic pressure
Optimal	Up to 120 mmHg	Up to 80 mmHg
Normal	Up to 130 mmHg	Up to 85 mmHg
High normal	130 – 139 mmHg	85 – 89 mmHg
Grade 1 hypertension	140 – 159 mmHg	90 – 99 mmHg
Grade 2 hypertension	160 – 179 mmHg	100 – 109 mmHg
Grade 3 hypertension	Over 180 mmHg	Over 110 mmHg

1. Introduction

Dear Customer, we are pleased that you have decided to purchase a blood pressure measurement device from HARTMANN. Tensoval duo control is a quality product for fully automatic blood pressure measurement on the human upper arm and is suitable for clinical and home application. Requiring no preliminary settings, this device conveniently inflates automatically for easy, quick and reliable measurement of the systolic and diastolic blood pressure as well as the pulse rate. In addition, you are informed if an irregular heart beat is detected.

We wish you all the best for your health.

2. General information on blood pressure

2.1 Significance of blood pressure values

To determine your blood pressure you need to measure two values:

- The systolic (upper) blood pressure: is produced when the heart contracts and pumps blood into the blood vessels.
- The diastolic (lower) blood pressure: This is the value measured when the heart muscle is dilated and again fills with blood.
- Blood pressure measurement values are indicated in mmHg.

The World Health Organisation (WHO) and the International Society for Hypertension (ISH) have developed the above mentioned classification for blood pressure values.

Established hypertension (high blood pressure) is defined as an office blood pressure reading of a systolic value of at least 140 mmHg and/or a diastolic value of at least 90 mmHg.

In general, low blood pressure (hypotension) is defined as a blood pressure value below 105 mmHg (systolic) and 60 mmHg (diastolic). However, this threshold between normal and too low blood pressure is not so accurately specified as the upward threshold in the direction of hypertension. You should also note that, unlike too high blood pressure values, too low blood pressure values are not usually expected to be associated with health risks.


2.2 Importance of self-measurement of blood pressure

Constantly elevated blood pressure multiplies the risk for other health problems. Physical consequences such as heart attack, stroke and organ damage are the most common causes of death worldwide. Daily blood pressure monitoring is thus an important measure which will help protect you from these risks.


3. Getting ready for self-measurement

3.1 Inserting / changing the batteries

Open the battery cover on the rear of the device by applying gentle pressure on the notch.

Insert the four batteries (type AA) so that the positive (+) and negative (-) polarities match the markings of the polarities on the battery compartment, as indicated. If the polarities do not match, the device will not work and the batteries may leak!

3.2 Setting date and time

Every time you insert batteries in the monitor, for the first time or for replacement, the monitor automatically switches to the date/time function. You can also re-set the date and time by pressing the START/STOP button for 5 seconds while the monitor is switched off.

When first used, the display shows the number "31" for the day and "12" for the month. The date set is thus the 31st December.

i You must correctly set the date/time before using the device for the first time, so that all memory and evaluation functions are working correctly.

The number on the left (date indication) in the display flashes. Press the \uparrow (+) or \downarrow (-) buttons to change the day displayed. For example, press \downarrow (-) twice to set the date to the 29th December. Store the current day by pressing the START/STOP button.

Now the number on the right flashes, which displays the month. You can set the month by pressing the \uparrow (+) or \downarrow (-) button and store by pressing the START/STOP button.

Now the year 2011 is displayed. You can change the year displayed as mentioned above by pressing the START/STOP button.

Then you may set the time. The number on the left in the display flashes, indicating 12:00 h. Once the desired number of hours has been set, store it by pressing the START/STOP button.

Now the number on the right flashes. Here you can change the minutes display and store it by pressing the START/STOP button.

i The measured values will remain in the memory when the batteries are changed. The date settings will also not be lost but the time settings must be reset.

3.3 10 golden rules for blood pressure measurement

1. Rest for approx. 5 minutes before measurement. Even deskwork increases blood pressure by an average of approx. 6 mmHg systolic and 5 mmHg diastolic.

2. Do not consume any nicotine or coffee up to one hour before measurement.

3. Do not measure when you have a strong urge to urinate. A full bladder can lead to an increase in blood pressure of approx. 10 mmHg.

4. Take measurements from the naked upper arm and while sitting comfortably upright.

5. In the case of using a wrist monitor, hold the cuff at heart level during the measuring procedure. The cuff of an upper arm monitor is automatically placed on the correct level on the arm.

6. Please do not talk or move during the measuring procedure. Talking increases the values by approx. 6 – 7 mmHg.

7. Wait at least one minute between two measurements, so that the vessels are relieved from pressure in preparation for a new measurement.

8. Enter the values in the blood pressure diary: note down the measured values, together with any drugs taken, date and time in your blood pressure diary.

9. Take measurements regularly. Even if your values have improved, you should continue to check them for monitoring purposes.

10. Always take measurements at the same time of day. Because a person has approx. 100,000 different blood pressure values every day, individual measurements have no significance. Only regular measurements at the same time each day over a longer period of time allow a meaningful assessment of blood pressure values.

Further information:

- You should take your blood pressure in a quiet place, in a relaxed and comfortable seated position.
- Measurement can be taken on the right or left arm. The arm giving higher readings should be used for long-term blood pressure monitoring.
- Always measure on the same arm and place your forearm in a relaxed position on a support.


3.4 Application of the cuff

Before you apply the cuff, insert the red cuff connector into the red cuff socket on the left side of the device.

Measurement must be taken on the naked upper arm. If your cuff is completely open, thread the end of the cuff through the metal hoop, so as to form a loop. In this instance, the Velcro fastener must be on the outside. Grasp the cuff by the grip tab (see illustration on the inside cover of the instructions for use) and wrap it over the upper arm.

4. Measuring the blood pressure

We recommend that you measure your blood pressure while sitting with your back supported by the back of the chair. Place both feet flat on the floor and the legs next to each other. Your legs should not be crossed. Place your forearm with the palm relaxed upwards on a support and make sure the cuff is at the level of the heart.



The recess (see illustration in the inside cover of the instructions for use) in the cuff, opposite the grip tab, should lie in the centre of the crook of the arm. The tube should lie in the crook of the arm, pointing towards the hand.

Now lightly bend your arm, take the free end of the cuff, wrap it firmly round your arm and close the Velcro fastener.

The cuff should be firm but not too tight. You should be able to push two fingers between the arm and the cuff. Make sure that the tube is not kinked or damaged.

Important: The correct application of the cuff is essential for the correct measurement result. The marking on the edge of the cuff helps you to choose the correct cuff size. The white arrow must point to an area within the sizing scale. If the white arrow is outside the scale, you need another cuff size (see Chapter 10 "Accessories and spare parts").

Do not switch the device on until the cuff has been applied, otherwise the cuff can become damaged through the resulting excess pressure. Press the START/STOP button.

The display of all display elements, followed by a flashing part of the progress bar, shows that the device is carrying out a self test and is ready for use.

Subsequently, automatic inflation begins. Three short beeps indicate that the necessary inflation pressure has been reached and that measurement is starting.

Then the measuring procedure begins with the release of air.

i If you require a higher inflation pressure, you can avoid having to repeat pumping by pressing the blue START/STOP button again after inflation starts and holding it down until the desired cuff pressure is reached. This should be approx. 30 mmHg over the systolic (upper) value.

You can track the measuring progress via the progress bar. This increases during the inflation phase and decreases in the measuring phase. During the measuring phase, you will also see the symbols for the two Duo Sensor Technology measuring methods. These indicate that the sensors are functioning correctly. In addition, the heart symbol shows your pulse beat, which is also measured.

 A long beep indicates the end of measurement. Then the systolic and diastolic blood pressure values appear simultaneously one below the other on the display as well as the pulse rate.

! Important: If you wish to stop the measuring procedure for any reason, simply press the START/STOP button. The inflation or measuring procedure is interrupted and an automatic pressure release occurs.

If you see this symbol  under the pulse rate, the device has detected an irregular heartbeat during the measurement. However, the measurement may also have been disrupted by body movement or speaking. It is best to repeat the measurement. If you regularly see this symbol during your blood pressure measurements, we recommend having your heart rhythm checked by your physician.

Once the measurement is finished  or  is displayed on the left.  represents the measured values for the first user. Using  you can store the measured values for a second user (see 5.1. Storing measured values). In order to switch off the device, press the START/STOP button. Otherwise the device will automatically switch off after 3 minutes.

5. Memory function

5.1 Storing measured values

The device features two memory buttons  and  to store measurement results for two different users.  represents the measured values for the first user,  represents the measured values for the second user. When a beep indicates the end of measurement, you may press  or  to assign the measured value to the respective person. This assignment can be made as long as the values are displayed. If you do not assign them, the measured value is automatically stored in the displayed measurement memory. The time the measurement was taken is always stored with the blood pressure values in order to determine e.g. the morning or evening averages correctly. The time stored in the device must therefore match the actual time of day (see chapter 3.2 “Setting date and time”).

5.2 Recalling the measured values

Tensoval duo control features the following measurement data storage (corresponding to the ESH - European Society of Hypertension guidelines)

- Single measurement memory
- Average value of all measured blood pressure values per user
- Average value of the morning values
- Average value of the evening values

The device must be switched off to recall data from the measurement memory. Press the memory button  for recalling the first user’s stored values and the memory button  for the second user’s values. The corresponding symbol  or  will be displayed. Tensoval duo control distinguishes between values measured in the morning and in the evening. Initially, the average of the morning values over the last seven days for the selected person is displayed  (“A” stands for average, 7 stands for seven days). After again pressing the memory button  or , the average of the evening values over the last seven days appears . After pressing the memory button again, the average of all measurements appears . The number of the stored measured values is displayed alternating with the sign “A”. Each time the  or  button is pressed after this, all stored measured values from the last 60 measurements appear

one after another, starting with the most recent value, along with the time, date and year, as well as displaying whether it is a morning or evening value   (01 is the latest measurement, 02 the next-to-last measurement, ...). The display alternates between the time, date and year of the measured value every 2 – 3 seconds.

i You can switch directly between the two users at any point in the memory display by pressing the respective memory button.


 Tensoval duo control stores per person ( or ) up to 60 measurements. The most recent measured value is always at memory position number 1. When all memory positions are occupied, the oldest value will be deleted.

i Important: Tensoval duo control follows the recommendations of the European Society of Hypertension (ESH) and distinguishes between the measured morning values and the measured evening values. This distinction is medically relevant, as the blood pressure changes in the course of the day. With this information, your physician is more likely to find the right treatment option for you in event of you receiving medical treatment for high blood pressure.

A measured value is assigned as a measured morning value if it is taken between 0:00 h and 12:00 h and as a measured evening value if it is taken between 12:01 h and 23:59 h. Please ensure that you always measure your blood pressure at the same time in the mornings and evenings.

The calculated general mean value is based on all measured values stored for the respective person. If only two measured values are in the memory, the mean value will be calculated from these two measured values. If only one measured value is in the memory, this measured value corresponds to the mean value.

 If there is no morning and/or evening value for the last seven days in the memory, then dashes are displayed instead of values in the display for morning and evening averages over the last seven days. If no values are stored at all in the memory, then dashes are also displayed instead of values in the display for the overall average.

i If an irregular heartbeat was detected during measurement, this information  is also stored and displayed when recalling the measured values from the device memory together with the systolic and diastolic blood pressure value, pulse rate, time, date and year.

You can cancel the memory recall at any time by pressing the START/STOP button. Otherwise, the device will switch off automatically after approx. 30 seconds. Even if the power supply fails, e.g. when changing the batteries, the stored values will still be available.

5.3 Deleting the measured values

You can delete the stored values for  and  separately. In order to delete all values of a user, press the memory button of the corresponding measurement memory  or . The morning average value for the last seven days will be displayed. If you press the button again and then hold down the memory button for 4 seconds, the numbers and signs except  or  start to flash in the display. If you hold down the memory button for another 4 seconds, all data for the selected user will be deleted. Now the display shows either  or .

If you would like to delete single values, please call up the corresponding single value to be deleted (see 5.2) and press the respective memory button for 4 seconds, so that the display flashes. After pressing for another 4 seconds, the respective single value is deleted. Now the display shows either  or .

⚠ Important: If you release the button ahead of time, no data will be deleted. When you delete an individual value, the next oldest measured value shifts to the memory position of the deleted measured value. When a measured value is deleted, the respective average values concerned are also recalculated.


5.4 Operation of guest mode

When a third person uses Tensoval duo control, it is recommended that they use the guest mode. This mode is used to store a measured value neither to the memory position  nor to .

This avoids bias to the series of measurements and the average values of the two main users of the device.

6. Explanation of error displays

Error which has occurred	Possible causes	Remedy
Device cannot be switched on	No batteries, batteries inserted incorrectly or batteries are dead.	Check batteries and insert four identical new batteries if necessary.
	Mains adapter not correctly connected or defective.	Ensure that the mains adapter is plugged into the connecting socket on the rear side of the device.
Cuff will not inflate	Cuff connector is incorrectly positioned in the connecting socket of the device.	Check connection between the red cuff connector and the red connecting socket.
	Wrong cuff type connected.	Check to make sure that only the approved Tensoval duo control cuffs and the corresponding connectors were used.
	The measuring signals could not, or not correctly, be read. This can be caused by incorrect application of the cuff, moving, talking or by a very weak pulse.	Check the correct position of the cuff. Do not talk or move during the measuring procedure. Also note the 10 golden rules and the information field from chapter 3.3

To take measurements in the guest mode, start the measuring procedure by simultaneously pressing the memory buttons  and . The START/STOP button need not be pressed. During and at the end of measurement, the display simultaneously shows the two symbols  and  next to the measured values. The measured value can thus not be applied to either person and the measured values will not be stored.

The device is also switched off in guest mode by pressing the START/STOP button. Otherwise the device will switch itself off automatically after 3 minutes.

Error which has occurred	Possible causes	Remedy
	Cuff does not inflate or does not inflate fast enough. This can, among other things, be due to a too loosely applied cuff or to movement during the measuring procedure.	Apply the cuff so that you are able to push about two fingers between the cuff and the upper arm. Air tube is not correctly inserted into the device. Check the correct position of the red connector plug. You should use a new cuff if this problem occurs often.
	Air discharge during the measuring procedure is too fast or too slow. The cuff could have become detached or loosened. Movement during the measuring procedure could also be the reason.	Check the correct position of the cuff. Do not move during the measuring procedure.
	Air tube is not correctly inserted into the device or kinked.	Check the correct position of the red connector plug or check the air tube for kinks.
	The pressure in the cuff exceeds 300 mmHg. Therefore, an automatic pressure release occurs.	Please rest for at least 1 minute and take the measurement again.
	If the battery symbol is flashing, the batteries are almost dead. Only a few more measurements are possible.	Keep a supply of batteries of the same type (type AA/LR06).
	If the battery symbol is permanently displayed, the batteries are dead and have to be replaced.	Insert new batteries of the same type (type AA/LR06). Then check the date / time and reset if necessary (see chapter 3.2).
Implausible measured values	Implausible measured values often occur due to inappropriate use of the device or mistakes during the measuring procedure.	Please note the 10 golden rules of blood pressure measurement and the safety information below. Then repeat the measurement.
	Wrong cuff size.	Use the cuff corresponding to your upper arm size.
	Cuff was placed on top of clothing.	Apply cuff on the naked skin.
	Rolled-up clothing impedes blood circulation.	Wear loose clothing. Make sure that rolled-up sleeves do not impede circulation in the upper arm.
	Cuff wrongly applied.	Take note of the instructions and images showing how to apply the cuff correctly to the upper arm.

Error which has occurred	Possible causes	Remedy
	Cuff tube folded or squashed.	Ensure that the cuff tube lies straight and loose.
	Cuff was not correctly inflated.	Check the correct position of the upper arm cuff.
	Moving, talking or excitement during the measuring procedure.	Please take measurements in a relaxed and comfortable position whilst seated. Do not talk or move during the measuring procedure.
	Lack of relaxation before taking a measurement.	Relax for at least 5 minutes before taking a measurement.
	Stimulants taken before the measurement.	Please avoid alcohol/ nicotine and caffeine for one hour before measurement.

Switch the device off if an error symbol appears. Check for all possible causes and note the 10 golden rules as well as the instructions for self-measurement from chapters 3 and 4. Relax for 1 minute and then repeat the measurement.

7. Explanation of symbols

 Consult instructions for use

 Caution, consult accompanying documents

 Protection from electric shock

 Manufacturer

 Disposal note

 Temperature limitation

 Humidity limitation

 Keep dry

8. Safety information concerning your health

8.1 Instructions for self-measurement

- Use the device exclusively for taking blood pressure measurements on the human upper arm.

- Do not take your blood pressure after taking a bath or exercising.
- Please wait one minute between two measurements.
- Single values are situation-dependent and have no significance.
- Even minor changes of inner and outer factors (e.g. deep breathing, stimulants, talking, excitement, climatic factors) lead to fluctuation in the blood pressure. This explains why deviating values are often measured by the physician or the pharmacist.
- The results of the measurement are basically dependent on the measuring site and position of the patient (sitting, standing, reclining). Furthermore, they are influenced by e.g. exertion and by the physiological preconditions of the patient.
- Do not leave the device unattended near toddlers or persons who cannot operate it themselves. There is a risk of strangulation through getting caught in the cuff tube. Swallowing small parts that have come off the device can also result in attacks of suffocation.
- Do not under any circumstances carry out blood pressure measurements on babies or toddlers.
- Please do not apply the cuff over a wound, as this may result in further injuries.

- If you have had a mastectomy, do not carry out the measurement on the arm on the affected side of the body.
- Please note that the pressure built up by the cuff can lead to temporary disruption to medical devices being simultaneously used on the same arm.
- If an intravenous treatment is being carried out or a venous catheter is present on the arm, blood pressure measurements can lead to injury. Please never use the cuff on the arm on which these conditions apply.
- If you are carrying out the measurement on another person, please ensure that the use of the blood pressure monitor does not result in persistent impairment of the blood circulation.
- Too frequent measuring within a short time period as well as continuous cuff pressure can interrupt the blood circulation and cause injury. Please make a pause between measurements and do not kink the air tube.

Consult your doctor before carrying out self-measurement of blood pressure if ...

- you are pregnant. Tensoval duo control has been clinically tested and can be used by pregnant women. In any case check with your physician whether and, if so, when you should carry out self-measurement, especially if you suffer from preeclampsia.
- you suffer from diabetes, hepatic disorders or narrowed blood vessels (for example, arteriosclerosis, peripheral arterial occlusive disease), because deviating measured values may occur in such cases.
- you suffer from certain blood diseases (e.g. haemophilia) or severely impaired blood flow, or if you take blood-thinning drugs.
- you suffer from severe cardiac arrhythmia.
- you wear a cardiac pacemaker. In this case, deviating measured values may occur. The blood pressure monitor itself has no impact on the cardiac pacemaker. Please note that the displayed pulse rate is not suitable for checking the rate of cardiac pacemakers.
- you tend to bruise and/or are sensitive to pain on pressure.

8.2 Arrhythmias, heart rhythm disorders

Heart rhythm disorders are disturbances of the normal rate or rhythm of the heartbeat. A distinction should be made between mild and severe heart rhythm disorders. This can only be determined by a special investigation carried out by a physician. Using Korotkoff technology, Tensoval duo control can measure correctly in the case of various types of heart rhythm disorders and thus provide correct measurement results.

 An irregular heartbeat is present when the heart rhythm deviates by more than 25 % from the average heart rhythm. The contraction of the heart muscle is stimulated by electrical signals. A disturbance in these electrical signals is called arrhythmia. This can be caused by genetic disposition, stress, age, lack of sleep, exhaustion, etc. A physician can establish whether an irregular heartbeat is the result of arrhythmia.

If this symbol  appears often, it can be a sign of heart rhythm disorders. In this case, consult your physician. Severe heart rhythm disorders may produce false measurements or impair the measuring accuracy. Please discuss with your physician whether self-measurement of blood pressure is suitable for you.

9. Maintenance of the device

Only use a soft, moistened cloth to clean the device. Please use neither a diluter nor alcohol, detergents or solvents. The cuff can be cleaned carefully with a lightly moistened cloth and mild soap solution. Do not completely immerse the cuff in water. Cleaning and disinfection of the cuff is recommended on a regular basis or after every use to prevent infection, especially in the case of application by several users. The disinfection, particularly on the inside of the cuff, should be performed by means of wipe disinfection. Use a disinfectant that is compatible with the cuff materials. To protect them from external influences, store the device and cuff as well as these instructions in the storage bag.

10. Accessories and spare parts

To ensure measuring accuracy, only use original HARTMANN accessories which may be obtained from your pharmacist or specialist medical supplier. Accurate measured results cannot be guaranteed outside the upper-arm circumference ranges specified here.

Upper arm circumference	Required cuff
17 – 22 cm	small
22 – 32 cm	medium
32 – 42 cm	large

Standard cuff, small, for upper arm circumference of 17 – 22 cm
Art. no. 900241

Standard cuff, medium, for upper arm circumference of 22 – 32 cm
Art. no. 900242

Standard cuff, large, for upper arm circumference of 32 – 42 cm
Art. no. 900243

Preformed cuff, medium, for upper arm circumference of 22 – 32 cm
Art. no. 900244

Tensoval mains adapter
Art. no. 900152

11. Transmission of blood pressure values to a computer

The product “Tensoval cardio control online”, which is available from medical retailers, enables the connection of the Tensoval duo control to the computer and the simple and secure transmission of the measured blood pressure values in an online application, meaning that you can

secure your measured values for a long period of time and read, store and print them from everywhere in graphic and tabular form.

An Internet connection is required for this purpose. Please note the hardware and software requirements if you purchase the Tensoval cardio control online.

Art. no. 900291, 900292, 900293

12. Warranty conditions

We give a 5-year warranty on this high-quality blood pressure monitor from the date of purchase and in accordance with the following conditions.

Claims must be made during the warranty period. The date of purchase has to be documented by the appropriately completed and stamped warranty certificate or proof of purchase.

Within the warranty period, HARTMANN shall replace or repair free-of-charge any device components which were caused to be faulty by material or manufacturing errors. This does not extend the warranty period.

The warranty does not cover damage resulting from improper use or tampering. The warranty will expire if the device is opened by an unauthorised person. Accessory parts that are subject to wear and tear (batteries, cuffs, mains adapters etc.) are excluded from the warranty. Claims for compensation are limited to the value of the merchandise; compensation for consequential damage or injury is expressly excluded.

In warranty cases, please send the device with cuff and, if applicable, the mains adapter together with the fully completed and stamped warranty certificate or proof of purchase direct, or via your dealer to the Customer Services department for your country.

13. Technical data

Type:	Automatic upper arm blood pressure monitor
Model:	Tensoval duo control
Measuring method:	Korotkoff and Oscillometric
Display range:	0 – 300 mmHg
Measuring range:	Systole (SYS): 50 – 250 mmHg Diastole (DIA): 40 – 160 mmHg Pulse: 40 – 160 pulse beats per minute
	The displaying of correct values outside the measuring range cannot be guaranteed.
Display unit:	1 mmHg
Technical measuring accuracy:	Cuff pressure: ± 3 mmHg, Pulse: $\pm 5\%$ of displayed pulse rate
Clinical measuring accuracy:	complies with the requirements of DIN EN 1060-4 and DIN EN ISO 81060-2; Korotkoff validation method: Phase I (SYS), Phase V (DIA)
Operating mode:	Continuous operation
Nominal voltage:	DC 6V
Power supply:	4 x 1.5 V Mignon alkaline-manganese (AA/LR06) batteries or optional HARTMANN Tensoval mains adapter
Battery capacity:	Tensoval duo control: > 1400 measurements Tensoval duo control Large: > 1400 measurements
Protection against electric shock:	Internally powered medical electrical (ME) equipment (when only batteries will be used). Applied part: Type BF  = Class II ME device (on application of the Tensoval mains adapter) IP20 (not protected against moisture)
Protection against harmful ingress of water or solid materials:	
Inflation pressure:	min. 140 mmHg
Comfort Air Technology:	Individual inflation pressure dependent on systolic blood pressure +30 mmHg
Automatic switch-off function:	3 minutes after end of measurement
Cuff:	Medium draw-clamp cuff 22 – 32 cm Large draw-clamp cuff 32 – 42 cm Small draw-clamp cuff 17 – 22 cm (optional) Medium preformed cuff 22 – 32 cm (optional)
Pressure release valve:	Electronically controlled linear valve
Memory capacity:	2 x 60 measurements and mean value over 7 days and overall mean value
Operating conditions:	Ambient temperature: +10 °C to +40 °C Relative humidity: 15 – 85 %, non-condensing Air pressure: 700 hPa – 1060 hPa
Storage/transport conditions:	Ambient temperature: –20 °C to +50 °C Relative humidity: 15 – 85 %, non-condensing
Serial number (SN):	inside the battery compartment
Computer interface to PC:	With the aid of Tensoval cardio control online, the stored measurement data can be readout and the measured values can be graphically displayed on your computer.

Reference to standards: DIN EN IEC 60601-1; DIN EN IEC 60601-1-2

PAUL HARTMANN AG
89522 Heidenheim, Germany

0 1 2 3

14. Power supply, disposal instructions

14.1 Batteries, mains adapter

- The four high-quality batteries included in the delivery guarantee at least 1400 measurements. Use only high-quality batteries (see specification in Chapter 13 "Technical data"). If you use less efficient batteries, we can no longer guarantee 1400 measurements.
- Never mix old and new batteries or batteries from different manufacturers.
- Remove exhausted batteries immediately.
- If you do not intend to use the device for a longer period of time, you should remove the batteries to prevent possible leakage.
- If you use the device with the mains adapter, position the device in such a way that you can disconnect the power supply at any time.

14.2 Disposal

- In the interests of environmental protection, exhausted batteries may not be disposed of in household waste. Please note the applicable waste disposal regulations or use public collecting bins.
-  This product is subject to European Directive 2012/19/EU on used electrical and electronic equipment and is marked accordingly. Never dispose of electronic equipment together with your household waste. Please obtain information concerning the local regulations on the proper disposal of electrical and electronic products. Proper disposal protects the environment and human health.

15. Safety information concerning the device

- This blood pressure monitor is not waterproof!
- This blood pressure monitor consists of high-quality electronic precision components. Protect the device from hard knocks or vibrations and do not drop it to the floor.

- Do not excessively bend or fold cuff and air tube.
- Never open the device. The device must not be altered, dismantled or repaired by the user. Repairs may only be carried out by authorised professionals.
- Never inflate the cuff when it is not properly applied to the upper arm.
- Use the device only with the approved upper arm cuff. Otherwise the device can be damaged externally or internally.
- The cuff tube may only be removed from the device by pulling on the red connector. Never pull on the tube itself!
- Do not expose the device to extreme temperatures or to humidity, dust or direct sunlight, as this may lead to malfunction.
- Keep the packaging, batteries and device out of reach of children.
- Please note the storage and operating conditions defined in Chapter 13 "Technical data". Storage or use outside the specified temperature and humidity ranges can have an impact on the measuring accuracy and functioning of the device.

Portable and mobile high-frequency and communication devices, such as telephone and mobile phone, can impair the functional capability of electronic medical devices and should therefore be kept at a distance of 3 m from the blood pressure monitor. In compliance with European Standard DIN EN IEC 60601-1-2 further information can be supplied by HARTMANN on request.

16. Legal requirements and guidelines

Tensoval duo control complies with the European regulations on which the provisions of the 93/42/EEC Medical Device Directive are based and bears the CE mark. The device complies, among others, with the requirements of the European Standard EN 1060: Non-invasive sphygmomanometers – Part 1:

General requirements and Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems, as well as with the standard EN 80601-2-30. The clinical testing of measuring accuracy was performed according to the European Standard EN 1060-4 as well as to EN 81060-2. The requirements of the ANSI/AAMI test protocol SP10-1992 are also fulfilled.

Over and beyond the legal requirements, the device has been clinically validated in accordance with the ESH-IP2 protocol of the European Society of Hypertension (ESH), the protocol of British Hypertension Society (BHS) and the protocol of the Deutsche Hochdruckliga (DHL).

17. Instructions for measurement function check

Each Tensoval duo control device has been thoroughly tested by HARTMANN for measuring accuracy and has been developed for a long useable service life. We recommend a measurement function check at intervals of 2 years in the case of **professionally used devices**, for example in pharmacies, medical practices or hospitals. You should also note the national regulations determined by the legislator. Measurement function checks can only be carried out either by competent authorities or authorized maintenance providers against compensation.

18. Instructions for the calibration mode

A device function check can be carried out on people or using a suitable simulator. Measurement function check involves testing for leak tightness of the pressure system and possible deviations of the pressure reading. Remove at least one battery to switch to calibration mode. Press down the START/STOP button and insert the battery again. Release the button after a few seconds and, after a few moments, two zeros will appear one above each other on the display. Instructions on the measurement function check will be supplied on request to competent authorities or authorised maintenance providers by HARTMANN.

19. Contact information for customer queries

AE – PAUL HARTMANN
Middle East FZE
P.O. Box 54525
Dubai
United Arab Emirates

AU – PAUL HARTMANN Pty. Ltd.
Level 6, 5 Ryder Boulevard
Rhodes, NSW 2138 Australia

CN – PAUL HARTMANN (Shanghai)
Trade Co., Ltd.
Shanghai 200233, P.R.China

GB – PAUL HARTMANN Ltd.
Heywood/Lancashire OL10 2TT

HK – PAUL HARTMANN
Asia-Pacific Ltd.
Hong Kong

ZA – HARTMANN South Africa
2194 Johannesburg

You will find more information on Tensoval duo control and other Tensoval products at www.tensoval.com

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■ English

- 1 Symbol for stored morning values
- 2 Symbol for stored evening values
- 3 Time / date
- 4 Systolic value
- 5 Diastolic value
- 6 Pulse
- 7 Irregular heartbeat
- 8 Battery symbol
- 9 Oscillometric measurement
- 10 Korotkoff's principle of blood pressure measurement
- 11 Memory user 2
- 12 Memory user 1
- 13 Progress bar